

AVSU COMPACT VERSION




The AVSU panel boards are conceived and realized according to the Directive 93/42/EEC and related ISO 7396 -1 norm.

AVSU groups are mainly used in low pressure medical gas networks as sectioning, control and supplying points, in case of maintenance or emergency.

The unit is available with pressure switch or pressure transducer according to the selected type of alarm.

Available up to 5 gas + vacuum.

AVSU STANDARD VERSION



The AVSU panel boards are conceived and realized according to the Directive 93/42/EEC and related ISO 7396 -1 norm.

AVSU groups are mainly used in low pressure medical gas networks as sectioning, control and supplying points, in case of maintenance or emergency.

The unit is available with pressure switch or pressure transducer according to the selected type of alarm.

Available up to 5 gas + vacuum.



STM SISTEMI TECNOLOGIE MEDICALI S.r.l. - VIA DELL'INFORMATICA, 7 - 20083 VIGANO DI GAGGIANO (MI) Italy
Tel: +39 02/90843021 - FAX: +39 02/90843022 - info@stmmedicali.eu - www.stmmmedicali.eu