


SINGLE STAGE CENTRAL STATION PANEL BOARD, type APR 180, WITH AUTOMATIC SWITCH-OVER DEVICE



The APR 180 CENTRAL STATION PANEL BOARD is conceived and manufactured according to MDD 93/42 EEC and subsequent amendments and harmonized ISO EN 7396-1, EN 10524-2, EN 737-3 standards.

Available for: O₂ / C.A. / N₂O / CO₂ / N₂

P. in: from 100 to 230 bar


P. out: 8 to 10 bar (factory pre settled at 8,5 bar)

Load: 180 m³/h


Built in signal board for ramp status

Varnished stainless steel support for a great increase of durability

Built in handles for easy installation and transport



EMERGENCY FIRST STAGE PANEL BOARD (THIRD SOURCE) 180 m³/h load



THE EMERGENCY FIRST STAGE PANEL BOARD (THIRD SOURCE) conceived and manufactured according with MDD 93/42 EEC and related ISO 7396-1 and UNI EN ISO 10524-2 norms.

Available for: O₂ / C.A. / N₂O / CO₂ / N₂

P. in: from 100 to 230 bar

P. out: 8 to 10 bar (factory pre settled at 8,5 bar)

Load: 180 m³/h

Varnished stainless steel support for a great increase of durability

