


## SINGLE STAGE CENTRAL STATION PANEL BOARD, type APR 40, WITH AUTOMATIC SWITCH-OVER DEVICE



**The APR 40 CENTRAL STATION PANEL BOARD is conceived and manufactured according to MDD 93/42 EEC and subsequent amendments and harmonized ISO EN 7396-1, EN 10524-2, EN 737-3 standards.**

Available for: O<sub>2</sub> / C.A. / N<sub>2</sub>O / CO<sub>2</sub> / N<sub>2</sub>

P. in: from 100 to 230 bar


P. out: 8 to 10 bar (factory pre settled at 8,5 bar)

Load: 40 m<sup>3</sup>/h


Built in signal board for ramp status

Varnished stainless steel support for a great increase of durability

Built in handles for easy installation and transport



## EMERGENCY FIRST STAGE PANEL BOARD (THIRD SOURCE) 40 m<sup>3</sup>/h load



**THE EMERGENCY FIRST STAGE PANEL BOARD (THIRD SOURCE)**  
conceived and manufactured according with MDD 93/42 EEC and  
related ISO 7396-1 and UNI EN ISO 10524-2 norms.

Available for: O<sub>2</sub> / C.A. / N<sub>2</sub>O / CO<sub>2</sub> / N<sub>2</sub>

P. in: from 100 to 230 bar

P. out: 8 to 10 bar (factory pre settled at 8,5 bar)

Load: 40 m<sup>3</sup>/h

Varnished stainless steel support for a great increase of durability

