

AFNOR MONOBLOC TERMINAL UNITS



AFNOR terminal units are conceived and realized according to the Directive 93/42/EEC and related ISO 7396-1, EN 9170-1, EN 737-1, EN 737-3 norms.

The terminal units can be mounted either in flush, surface box as well as in bed-head units and ceiling pendants.

They are composed by:


- Nickel-plated gas specific terminal unit realized according to the EN 9170 and AFNOR NFS 90-116 norms



GAS	VACUUM
A	Ø 7
B	Ø 8,1
GROOVES N°	2



GAS	NITROUS OXIDE
A	Ø 7
B	Ø 7,1
GROOVES N°	4



GAS	OXYGEN
A	Ø 7
B	Ø 8,1
GROOVES N°	3





GAS	MEDICAL AIR
A	Ø 8
B	Ø 7,1
GROOVES N°	2



GAS	CARBON DIOXIDE
A	Ø 5
B	Ø 7,1
GROOVES N°	5



GAS	NITROGEN
A	Ø 5,5
B	Ø 7,1
GROOVES N°	4



SURFACE MOUNTING BOX



MAINTENANCE KIT Kit for maintenance all in one

