

FLOWMETER FOR OXYGEN AND MEDICAL AIR



The range of flowmeters manufactured by STM is conceived and realized according to MDD 93/42/EEC and related norms.

Flowmeters are available both single and double solution.

Scale: 0-6, 10, 15, 30 lt/mnt

Probes: UNI/AFNOR/DIN/BS

Downstream flow regulation and floater for extreme precise measuring.



HUMIDIFIER



300 cc humidifier composed by:

Nickel-plated brass cap with tube holder nipple and safety valve,
stainless steel diffuser, makrolon jar sterilizable at 121°C for 15 minutes
Weight Kg 0,610



STM SISTEMI TECNOLOGIE MEDICALI S.r.l. - VIA DELL'INFORMATICA, 7 - 20083 VIGANO DI GAGGIANO (MI) Italy
Tel: +39 02/90843021 - FAX: +39 02/90843022 - info@stmmedicali.eu - www.stmmedicali.eu