

**DOUBLE STAGE CENTRAL STATION PANEL BOARD,  
type APR 100 HD, WITH AUTOMATIC SWITCH-OVER DEVICE  
BY MEANS OF SOLENOID VALVES (N.O.)**



**The APR 100 HD CENTRAL STATION PANEL BOARD is conceived and manufactured according to MDD 93/42 EEC and subsequent amendments and harmonized ISO EN 7396-1, EN 10524-2, EN 737-3 standards and designed to meet the requirements of HTM 02-01**

Available for: O<sub>2</sub> / C.A. / N<sub>2</sub>O / CO<sub>2</sub> / N<sub>2</sub>

P. in: from 100 to 230 bar

P. out: 4 to 8 bar (factory pre settled at 5 bar)


Load: 100 m<sup>3</sup>/h

Built in signal board for ramp status

Varnished stainless steel support for a great increase of durability

Built in handles for easy installation and transport

Fully computerized management system with 7" LCD touch screen



## ERM



**THE ERM is conceived and manufactured according to MDD 93/42 EEC and subsequent amendments and harmonized ISO EN 7396-1, EN 10524-2, EN 737-3 standards and designed to meet the requirements of HTM 02-01**

Available for: O<sub>2</sub> / C.A. / N<sub>2</sub>O / CO<sub>2</sub> / N<sub>2</sub>

P. in: from 100 to 230 bar

P. out: 4 to 8 bar (factory pre settled at 5 bar)

Load: 100 m<sup>3</sup>/h

Varnished stainless steel support for a great increase of durability

