


MULTIPLE PRESSURE REDUCER PANELS, WR MODEL



The WR model pressure reducer and multiple panel boards, are conceived and realized according to MDD 93/42/EEC and related ISO 7396-1 and ISO 10524-2 norms.

The multiple pressure reducer panel boards are conceived in order to satisfy the specific demands of each area of a medical gas network having the possibility to carry up to 5 GAS+VACUUM as standard configuration.

MULTIPLE PRESSURE REDUCER PANELS, WR MODEL WITH INTEGRATED ALARM



MULTIPLE PRESSURE REDUCER PANELS, WR MODEL WITHOUT ALARM

