


**PRESSURE REDUCER WITH CLICK ADJUSTMENT FLOW**

**stm**  
SISTEMI TECNOLOGIE MEDICALI S.P.A.

- Pressure reducer for cylinder, preset at 3,5 bar
- Chrome – brass body
- Teflon membrane
- Safety valve
- Inlet fitting UNI – NF – DIN – BS
- Flow rate 0 - 15 Lt/min
- Intervals 1-1,5-2-3-4-6-9-12-15 Lt/min
- High pressure gauge 0-315 / 0-160 bar
- Output fitting 9/16 UNF with hose connection D.6
- Max load 15 cu.mt/hr
- Temperature range – 20 ° + 50 ° C
- Weight Kg. 0,690

## PRESSURE REDUCER WITH CLICK ADJUSTMENT FLOW

**PRESSURE REDUCER WITH CLICK ADJUSTMENT FLOW**

Pressure reducer for cylinder, preset at 3,5 bar

Chrome – brass body

Teflon membrane

Safety valve

Inlet fitting UNI – NF – DIN – BS

Flow rate 0 - 15 Lt/min

Intervals 1-1,5-2-3-4-6-9-12-15 Lt/min

High pressure gauge 0-315 / 0-160 bar

Output fitting 9/16 UNF with hose connection D.6

Max load 15 cu.mt/hr

Temperature range – 20 ° + 50 ° C

Weight Kg. 0,630



**FLOWMETER** with snap adjustment flow regulation, for Oxygen or Medical Air, working pressure 3,5 bar

Chrome – brass body

Inlet fitting UNI – NF – DIN – BS

Flow rate 0 - 15 Lt/min

Intervals 1-1,5-2-3-4-6-9-12-15 Lt/min

On request

Flow rate 0 a 6 Lt/min

Intervals 0,5-1-1,5-2-2,5-3-4-5-6

Gauge 0-10 bar

Output fitting 9/16 UNF with hose connection D.6

Temperature range – 20 ° + 50 ° C

Weight Kg. 0,780